

**ROTANODE™  
XRR-4631G**

**Rotating Anode X-ray Tube Assembly**

- ◆ High speed rotating anode X-ray tube assembly for high energy radiographic operations.
- ◆ For the purpose of general diagnostic X-ray procedures.
- ◆ This tube has foci 1.2 and 0.6, and is available for a maximum tube voltage 150kV.
- ◆ This tube assembly has specially processed rhenium-tungsten faced molybdenum target of 100mm diameter anode disc and is accommodated with IEC60526 high-voltage cable receptacles.

**General Data**

**IEC Classification (IEC60601-1:2005+A1:2012) ..... Class I ME EQUIPMENT**

**Electrical:**

Circuit:

High Voltage Generator ..... Constant Potential High-voltage Generator  
 Grounding ..... Center-grounded

Nominal X-ray Tube Voltage:

Radiographic ..... 150 kV

Nominal Focal Spot Value:

Large Focus ..... 1.2

Small Focus ..... 0.6

Nominal Anode Input Power (at 0.1s):

Large Focus ..... 180 Hz    60 Hz

Small Focus ..... 100 kW    58 kW

Large Focus ..... 40 kW    23 kW

Nominal Radiographic Anode Input Power:

Large Focus ..... 180 Hz    60 Hz

Small Focus ..... 94 kW    56 kW

Large Focus ..... 36 kW    23 kW

★The information contained herein is presented only as a guide for the application of our products. No responsibility is assumed by Canon Electron Tubes & Devices Co., Ltd. (CETD) for any infringements of patents or other rights of the third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of CETD or others.

★The information contained herein may be changed without prior notice. It is therefore, advisable to contact to CETD before processing with the design of equipment incorporating this product.

Motor Ratings:

Stator: XS-BM

	Starting		Running	
	180	60	180	60
Driven Frequency [Hz]	180	60	180	60
Input Power [W]	3600	1760	320	240
Voltage <sup>3) 7)</sup> [V]	260	140	70	50
Current <sup>4)</sup> [A]	15	13.5	5	5
Min. Speed Up <sup>1) 7)</sup> [s]	1.3	1.3	-	-
Capacitor [ $\mu$ F]	20	65	20	65
Min. Braking <sup>2) 7)</sup> [s]	3.0 / 60V (DC)			

- Note 1) The speed up time from normal speed to high speed is 2/3 times of the specified speed up time from 0 to high speed, which is described on motor rating table.  
 2) To be applied for high speed rotation.  
 3) Applied voltage between common and main terminal.  
 4) Common current.  
 5) The every applied voltage must be never exceeded 110% of the above specification.  
 6) No more than two high speed starts per minute are permissible.  
 7) The speed-up time is allowed up to 110% of the above specification.  
 8) The generator manufacture may choose different values. The above table is one of the recommend conditions.

Anode Speed\*:

180 Hz ..... Minimum 9700 min<sup>-1</sup>  
 60 Hz ..... Minimum 3120 min<sup>-1</sup>

Stator Resistance:

Common-Main Winding ..... 6  $\Omega$   
 Common-Auxiliary Winding ..... 11  $\Omega$

Resistance between Housing and Low Voltage Terminals ..... Minimum 2 M $\Omega$

Normal Operating Range of the Housing Temperature ..... 16 ~ 75 °C

Mode of Operation ..... Intermittent

**Mechanical:**

Dimensions .....	See Dimensional Outline
Overall Length .....	445 mm
Maximum Diameter .....	166 mm
Target:	
Anode Angle .....	12 degrees
Diameter .....	100 mm
Construction .....	Rhenium-Tungsten-faced Molybdenum
Filtration:	
Permanent Filtration .....	1.1 mm AL / 75 kV IEC60522:1999
Available Additional Filter combination (0.4 - 1.5 mm) .....	Maximum 2.6 mm AL / 75 kV
Radiation Protection (In accordance with IEC60601-1-3:2008):	
Leakage Technique Factor .....	150 kV, 3.4 mA
X-ray Coverage .....	430 × 430 mm at SID 1000 mm
Weight (Approx.) .....	20 kg
High Voltage Receptacle .....	To meet the requirements of IEC60526 Corrigendum1:2010
Cooling Method .....	Natural or Forced Air
Tube Housing Model Number .....	XH-1019

## Absolute Maximum and Minimum Ratings

(At any time, these values must not be exceeded.)

### Maximum X-ray Tube Voltage:

Radiographic .....	150 kV
Between Anode (or Cathode) and Ground .....	75 kV
Minimum X-ray Tube Voltage .....	40 kV
Maximum X-ray Tube Current .....	See Rating Charts
Large Focus .....	1000 mA
Small Focus .....	400 mA

### Maximum Filament Current:

Large Focus .....	5.2 A
Small Focus .....	5.1 A

### Filament Voltage:

Large Focus (At Maximum Filament Current 5.2 A) .....	11.3 ~ 15.4 V
Small Focus (At Maximum Filament Current 5.1 A) .....	5.9 ~ 8.1 V

Filament Frequency Limits ..... 0 ~ 25 kHz

Continuous Anode Input Power ..... 300 W (420 HU/s)

### Thermal Characteristics:

Anode Heat Content .....	285 kJ (400 kHU)
Maximum Anode Heat Dissipation .....	1200 W (1690 HU/s)
X-ray Tube Assembly Heat Content .....	1480 kJ (2000 kHU)

### Nominal Continuous Input Power:

Without Air-circulator .....	222 W (18 kHU/min)
------------------------------	--------------------

## Environmental Limits

### Operating Limits:

Temperature .....	10 ~ 40 °C
Humidity .....	30 ~ 85 %
	(No Condensation)
Atmospheric Pressure .....	70 ~ 106 kPa


### Shipping and Storage Limits:

Temperature .....	-20 ~ 70 °C
Humidity .....	20 ~ 90 %
	(No Condensation)
Atmospheric Pressure .....	50 ~ 106 kPa


## Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage  
Constant Potential High-voltage Generator  
Stator Power Frequency 180Hz

Nominal Focal Spot Value: 1.2 ■




Nominal Focal Spot Value: 0.6 □




## Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage  
Constant Potential High-voltage Generator  
Stator Power Frequency 60Hz

Nominal Focal Spot Value: 1.2 ■




Nominal Focal Spot Value: 0.6 □



## Emission & Filament Characteristics


Constant Potential High-Voltage Generator

Nominal Focal Spot Value: 1.2 ■



For Reference Only


Nominal Focal Spot Value: 0.6 □




For Reference Only

## Thermal Characteristics

X-ray Tube Assembly Heating / Cooling Curve



Anode Heating / Cooling Curve




The heating curves are showing examples of average input power to the anode in operation.






## Dimensional Outline

Unit mm



Note) Do not connect terminal No. 1 and No. 5 or No. 6 in series circuit.



### EXPLANATION OF SYMBOLS

CATHODE TERMINAL

C : COMMON

L : LARGE FOCUS

S : SMALL FOCUS

### TERMINAL CONNECTIONS


C1 : COMMON

M : MAIN WINDING OF THE STATOR

A : AUX. WINDING OF THE STATOR

NC : NON-CONNECTION

ET : EARTH TERMINAL



▲ : CENTRAL X-RAY  
ANODE & CATHODE TERMINAL  
: IEC60526 TYPE



**CANON ELECTRON TUBES & DEVICES CO., LTD.**

Marketing Engineering Group, Sales Department  
1385, Shimoishigami, Otawara-shi, Tochigi 324-8550, Japan  
Tel: +81-287-26-6666 Fax: +81-287-26-6060  
<https://etd.canon>

·The head office of Canon Electron Tubes & Devices Co., Ltd. has been certified to meet all the requirements of Environmental Management System ISO14001.  
·Canon Electron Tubes & Devices Co., Ltd. has been certified to meet all the requirements of Quality Management Systems ISO9001 and ISO13485.  
Product scope is referred to the following URL. <https://etd.canon/eng/company/quality.htm>.