

**ROTANODE™  
XRR-6652X**

**Rotating Anode X-ray Tube Assembly**

- ◆ High speed rotating anode X-ray tube assembly for high energy radiographic operations.
- ◆ For the purpose of general diagnostic X-ray procedures.
- ◆ This tube has foci 0.8 and 0.3, and is available for a maximum tube voltage 150 kV.
- ◆ The anode disk is constructed with specially processed rhenium-tungsten-faced molybdenum target of 100 mm diameter.

**General Data**

**IEC Classification (IEC60601-1:2005+A1:2012) ..... Class I ME EQUIPMENT**

**Electrical:**

Circuit:

High Voltage Generator ..... Constant Potential High-Voltage Generator  
 Grounding ..... Center-Grounded

Nominal X-ray Tube Voltage:

Radiographic ..... 150 kV  
 Fluoroscopic ..... 125 kV

Nominal Focal Spot Value:

Large Focus ..... 0.8  
 Small Focus ..... 0.3

Nominal Anode Input Power (at 0.1s) ..... See rating charts

	180 Hz	60 Hz	50 Hz
Large Focus .....	52 kW	30 kW	29 kW
Small Focus .....	12 kW	6.5 kW	6.1 kW

Nominal Radiographic Anode Input Power:

	180 Hz	60 Hz	50 Hz
Large Focus .....	52 kW	30 kW	29 kW
Small Focus .....	12 kW	6.5 kW	6.1 kW

★The information contained herein is presented only as a guide for the application of our products. No responsibility is assumed by Canon Electron Tubes & Devices Co., Ltd. (CETD) for any infringements of patents or other rights of the third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of CETD or others.  
 ★The information contained herein may be changed without prior notice. It is therefore, advisable to contact to CETD before processing with the design of equipment incorporating this product.

Motor Ratings:

Stator: XS-AG

	Starting		Running	
	180	60	180	50/60
Driven Frequency [Hz]	180	60	180	50/60
Input Power [W]	3710	800	200	90
Voltage <sup>3) 5)</sup> [V]	420	190	80	48
Current <sup>4)</sup> [A]	9.6	10.2	3.0	2.2
Min. Speed Up <sup>1) 7)</sup> [s]	1.2	0.8	-	-
Capacitor [ $\mu$ F]	6	44	6	44
Min. Braking <sup>2) 7)</sup> [s]	2.0 (DC 100V)			

- Note 1) The speed up time from normal speed to high speed is 2/3 times of the specified speed up time from 0 to high speed, which is described on motor rating table.  
 2) To be applied for high speed rotation.  
 3) Applied voltage between common and main terminal.  
 4) Common current.  
 5) The every applied voltage must be never exceeded 110% of the above specification.  
 6) No more than two high speed starts per minute are permissible.  
 7) The speed-up time is allowed up to 110% of the above specification.

Anode Speed:

- 180 Hz ..... Minimum 9700 min<sup>-1</sup>  
 60 Hz ..... Minimum 3200 min<sup>-1</sup>  
 50 Hz ..... Minimum 2700 min<sup>-1</sup>

Stator Resistance:

- Common-Main Winding ..... 9.4  $\Omega$   
 Common-Auxiliary Winding ..... 28.3  $\Omega$   
 Resistance between Housing and Low Voltage Terminals ..... Minimum 2 M $\Omega$   
 Normal Operating Range of the Housing Temperature ..... 16 ~ 75 °C  
 Mode of Operation ..... Intermittent

**Mechanical:**

Dimensions .....	See dimensional outline
Overall Length .....	445.3 mm
Maximum Diameter .....	165.6 mm
Target:	
Anode Angle .....	12 degrees
Diameter .....	100 mm
Construction .....	Rhenium-Tungsten faced Molybdenum
Permanent Filtration .....	1.5 mm Al / 75 kV IEC60522:1999
Radiation Protection (In accordance with IEC60601-1-3:2008):	
Leakage Technique Factor .....	150 kV, 3.3 mA
X-ray Coverage .....	430 × 430 mm at SID 1000 mm
Weight:	
Tube Housing Unit .....	Approx. 23 kg
Heat Exchanger, Hose .....	Approx. 17 kg
High Voltage Receptacle .....	To meet the requirements of IEC60526 Corrigendum1:2010
Cooling Method .....	Oil Circulation
Tube Housing Model Number .....	XH-1022

## **Absolute Maximum and Minimum Ratings**

**(At any time, these values must not be exceeded.)**

**Maximum X-ray Tube Voltage:**

Radiographic .....	150 kV
Fluoroscopic .....	125 kV

Between Anode (or Cathode) and Ground ..... 75 kV

Minimum X-ray Tube Voltage ..... 40 kV

**Maximum X-ray Tube Current:**

Large Focus .....	700 mA
Small Focus .....	150 mA

**Maximum Filament Current:**

Large Focus .....	5.5 A
Small Focus .....	5.0 A

**Filament Voltage:**

Large Focus (At maximum filament current 5.5 A) .....	11.6 ~ 14.8 V
Small Focus (At maximum filament current 5.0 A) .....	4.6 ~ 5.8 V

Filament Frequency Limits ..... 0 ~ 25 kHz

Continuous Anode Input Power ..... 500 W (700 HU/s)  
(Fluoroscopic, Radiographic or mixed exposure)

**Thermal Characteristics:**

Anode Heat Content .....	420 kJ (600 kHU)
Maximum Anode Heat Dissipation .....	1670 W (2300 HU/s)
X-ray Tube Assembly Heat Content .....	1420 kJ (2000 kHU)

**Nominal Continuous Input Power:**

Oil Circulation (Heat Exchanger, Without Tube Cover) ..... 1000 W (84.7 kHU/min)

## Environmental Limits

### Operating Limits:

Temperature .....	10 ~ 40 °C
Humidity .....	30 ~ 85 %
	(No condensation)
Atmospheric Pressure .....	70 ~ 106 kPa


### Shipping and Storage Limits:

Temperature .....	-20 ~ 70 °C
Humidity .....	20 ~ 90 %
	(No condensation)
Atmospheric Pressure .....	50 ~ 106 kPa


## Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 180Hz

Nominal Focal Spot Value: 0.8 ■




Nominal Focal Spot Value: 0.3 □




## Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 60Hz

Nominal Focal Spot Value: 0.8 ■




Nominal Focal Spot Value: 0.3 □




## Maximum Rating Charts (Absolute Maximum Rating Charts)

Conditions: Tube Voltage  
Constant potential high-voltage generator  
Stator Power Frequency 50Hz

Nominal Focal Spot Value: 0.8 ■



Nominal Focal Spot Value: 0.3 □






## Emission & Filament Characteristics


Constant potential high-voltage generator

Nominal Focal Spot Value: 0.8 ■



For Reference Only


Nominal Focal Spot Value: 0.3 □




For Reference Only

## Thermal Characteristics

X-ray Tube Assembly Heating / Cooling Curve




Anode Heating / Cooling Curve



### Dimensional Outline

Unit: mm



### Dimensional Outline (Heat Exchanger)

Unit: mm

