

**X-RAY TUBE**  
**D-125**  
**D-125S**  
**D-125SB**

**Stationary Anode X-ray Tube**

- ◆ Especially designed for dental X-ray unit.
- ◆ Provided with on insulation cylinder (D-125S) and lead cylinder (D-125SB).
- ◆ These tubes have focus 1.2, and are available for maximum tube voltage 100 kV or 125 kV.
- ◆ Installed in the same enclosure with the high voltage transformer.

**General Data**

**Electrical:**

Circuit:

High Voltage Generator ..... Constant Potential High-Voltage Generator  
 Grounding ..... Center-Grounded

Nominal X-ray Tube Voltage:

D-125 ..... 125 kV  
 D-125S, D-125SB ..... 100 kV

Nominal Focal Spot Value ..... 1.2

Nominal Anode Input Power (at 1.0s) (See rating charts) ..... 2700 W

Nominal Radiographic Anode Input Power ..... 2900 W

Exposure Duty Cycle ..... 1:60 or more  
 (Exposure Time : Interval Time)

**Mechanical:**

Dimensions:

Overall Length ..... See dimensional outline  
 Maximum Diameter ..... See dimensional outline

Target:

Anode Angle ..... 16 degrees  
 Material ..... Tungsten

Inherent Filtration ..... At least 0.8 mm Al at 50 kV

X-ray Coverage ..... 354 × 354 mm at SID 750 mm

★The information contained herein is presented only as a guide for the application of our products. No responsibility is assumed by Canon Electron Tubes & Devices Co., Ltd. (CETD) for any infringements of patents or other rights of the third parties which may result from its use.  
 No license is granted by implication or otherwise under any patent or patent rights of CETD or others.  
 ★The information contained herein may be changed without prior notice. It is therefore, advisable to contact to CETD before processing with the design of equipment incorporating this product.

Weight (approx.):

D-125 .....	350 g
D-125S .....	430 g
D-125SB .....	780 g

Cooling Method ..... Oil immersed (60°C Max.) and convection oil cooling.

Tube Holding:

D-125 .....	Holding the glass envelope of the anode end and cathode end, or the screw of the anode shank.
D-125S .....	Holding the insulation cylinder.
D-125SB .....	Holding the Insulation cylinder or the lead cylinder.

**Absolute Maximum and Minimum Ratings**  
(At any time, these values must not be exceeded.)

Maximum X-ray Tube Voltage:

D-125 .....	125 kV
D-125S, D-125SB .....	100 kV

Between Anode (or Cathode) and Ground:

D-125 .....	62.5 kV
D-125S, D-125SB .....	50 kV

Minimum X-ray Tube Voltage ..... 40 kV

Maximum X-ray Tube Current ..... 40 mA

Maximum Filament Current ..... 3.1 A

Filament Voltage:

At maximum filament current (3.1A) ..... 2.7 ~ 3.6 V

Filament Frequency Limits ..... DC or AC (Sine Wave) 0 ~ 20 kHz

Thermal Characteristics:

Anode Heat Content ..... 35 kJ

Maximum Anode Heat Dissipation ..... 250 W

Maximum Radiographic Exposure Time ..... 10 s

**Environmental Limits**

Operating Limits (in dielectric oil):

Oil Temperature .....	10 ~ 60°C
Oil pressure .....	70 ~ 140 kPa

Shipping and Storage Limits:

Temperature .....	-40 ~ 70°C
Humidity .....	10 ~ 90 %


(No condensation)

Atmospheric Pressure ..... 50 ~ 106 kPa

## Maximum Rating Charts (Absolute maximum rating charts)


Constant Potential High-Voltage Generator

Nominal Focal Spot Value: 1.2




## Emission & Filament Characteristics

Constant Potential High-Voltage Generator




This graph indicates typical characteristics.

### Anode Heating / Cooling Curve




### Dimensional Outline of D-125

Unit: mm





### Dimensional Outline of D-125S

Unit: mm



### Dimensional Outline of D-125SB

Unit: mm



F : FILAMENT  
 C1, C2 : CATHODE  
 N : NOT USE



**CANON ELECTRON TUBES & DEVICES CO., LTD.**

Marketing Engineering Group, Sales Department  
1385, Shimoishigami, Otawara-shi, Tochigi 324-8550, Japan  
Tel: +81-287-26-6666 Fax: +81-287-26-6060  
<https://etd.canon>

·The head office of Canon Electron Tubes & Devices Co., Ltd. has been certified to meet all the requirements of Environmental Management System ISO14001.  
·Canon Electron Tubes & Devices Co., Ltd. has been certified to meet all the requirements of Quality Management Systems ISO9001 and ISO13485.  
Product scope is referred to the following URL. <https://etd.canon/eng/company/quality.htm>.