

**ANALIX™  
E7340X**

- ◆ The ANALIX E7340X, a water cooled, grounded cathode X-ray tube designed for X-ray fluorescence spectrometry.
- ◆ The maximum input power 3.0 kW.
- ◆ Suitable for the microanalysis of light metals such as Ca and K contained in wasted water.

**General Data**

**Electrical:**

Circuit:	
High Voltage Generator .....	Constant Potential High-Voltage Generator
Grounding .....	Cathode Grounded
Nominal X-ray Tube Voltage .....	60 kV
Focal Spot .....	ϕ 14 mm
Starting Tube Voltage .....	20 kV or Less

**Mechanical:**

Dimensions:	
Overall Length .....	459 mm
Diameter of Maximum Part .....	ϕ 124 mm
Target:	
Angle .....	90 Degrees to The Tube Axis.
Material .....	Rh
Inherent Filtration .....	127 µm Beryllium
Weight .....	Approx. 9.5 kg
Cooling:	
Anode Target .....	Water cooled
Water Quantity .....	Minimum 2.4 ℓ/min
Water .....	Demineralized, Deoxygenated and High Resistibility Water
Water Resistibility .....	Minimum 500 kΩ·cm
Maximum Inlet Temperature .....	45°C
Maximum Water Pressure .....	490 kPa
Pressure Drop .....	68.6 kPa at 2.4 ℓ/min
Cathode Jacket .....	Water Cooled
Water Quantity .....	Minimum 2.0 ℓ/min
Water .....	Tap Water
Maximum Inlet Temperature .....	55 °C
Maximum Water Pressure .....	490 kPa
Pressure Drop .....	49 kPa at 2.0 ℓ/min

★The information contained herein is presented only as a guide for the application of our products. No responsibility is assumed by Canon Electron Tubes & Devices Co., Ltd. (CETD) for any infringements of patents or other rights of the third parties which may result from its use.  
No license is granted by implication or otherwise under any patent or patent rights of CETD or others.  
★The information contained herein may be changed without prior notice. It is therefore, advisable to contact to CETD before processing with the design of equipment incorporating this product.

## Absolute Maximum and Minimum Ratings

(At any time, these values must not be exceeded.)

Maximum X-ray Tube Voltage .....	60 kV
Minimum X-ray Tube Voltage .....	20 kV
Maximum X-ray Tube Current .....	100 mA
Maximum Filament Current .....	10 A
Filament Frequency Limits .....	20 kHz or Less
Maximum Input .....	3.0 kW

## Environmental Limits

### Operating Limits:


Temperature .....	10 ~ 40 °C
Humidity .....	85 % or Less (No condensation)
Atmospheric Pressure .....	70 ~ 106 kPa
Vibration and Shock .....	Use it in environment without vibration and a shock

### Shipping and Storage Limits:

Temperature .....	-20 ~ 70 °C
Humidity .....	85 % or Less (No condensation)
Atmospheric Pressure .....	50 ~ 106 kPa
Vibration .....	4.9 m/s <sup>2</sup> or Less
Shock .....	98 m/s <sup>2</sup> or Less


## Maximum Rating Charts (Absolute Maximum Rating Charts)

Constant Potential High-Voltage Generator



## Emission & Filament Characteristics


Constant Potential High-Voltage Generator



Note: This graph indicates typical characteristics.

### Dimensional Outline (1)

Unit: mm




- Note:
1. Mark  $\uparrow$  indicates the direction of central X-ray.
  2. Filamet cable terminal: Amp power lock connector.  
(No. 108-11026, UL No. E28476) Made by AMP (Japan)
  3. High voltage cable receptacle: IEC60526 Type.

### Dimensional Outline (2)

Unit: mm

Ⓐ Enlargement



---

 **CANON ELECTRON TUBES & DEVICES CO., LTD.**  
Marketing Engineering Group, Sales Department  
1385, Shimoishigami, Otawara-shi, Tochigi 324-8550, Japan  
Tel: +81-287-26-6666 Fax: +81-287-26-6060  
<https://etd.canon>

·The head office of Canon Electron Tubes & Devices Co., Ltd. has been certified to meet all the requirements of Environmental Management System ISO14001.  
·Canon Electron Tubes & Devices Co., Ltd. has been certified to meet all the requirements of Quality Management Systems ISO9001 and ISO13485.  
Product scope is referred to the following URL. <https://etd.canon/eng/company/quality.htm>.